

## A cell history recorder for the incubator

**+ Compact design and stability: ideal for use in incubators - even several devices at the same time**

**+ Image processing algorithms for cell counting and well coverage**

**+ 24 digital microscopes with autofocus for sharp images - various adapter for different cell culture vessels**

**+ USB 3.0 interface to computer**

## The OWL for your laboratory

- ▶ High Level of perception like the OWL in the nature.
- ▶ 24-channel microscope for rapid and **automated cell culture microscopy**.
- ▶ **Sharp images** of your cells even with incubator use.



## For which application?

Documentation and determination of **cell growth**.◀

Determine **influences** on your cell culture by compound addition.◀

**Migration** assays.◀

Stem cell observation.◀



## Why choose the OWL?

- ▶ Integrated image processing algorithms to evaluate your **long-term experiment**.
- ▶ Provides **fast and accurate information** about the current state of the sample.
- ▶ **No cell stress** because of no changing condition/ temperature.
- ▶ **Reliable documentation** of your experiment.



[www.labtastic.shop](http://www.labtastic.shop)

A product of innoME GmbH - In der Tütenbeke 36 - 32339 Espelkamp - 05772/ 560 500